

Response to Article “Oriented Graphene Oxide Scaffold Promotes Nerve Regeneration in vitro and in vivo” [Letter]

Novrynda Eko Satriawan^{1,*}, Ariyani Noviantari^{1,2,*}

¹Doctoral Program in Biomedical Science, Faculty of Medicine, Universitas Indonesia, Jakarta, Indonesia; ²Center for Biomedical Research, Research Organization for Health, National Research and Innovation Agency (BRIN), Cibinong Science Center, Cibinong-Bogor, West Java, Indonesia

*These authors contributed equally to this work

Correspondence: Novrynda Eko Satriawan, Doctoral Program in Biomedical Science, Faculty of Medicine, Universitas Indonesia, Jalan Salemba Raya No. 6, Jakarta, Indonesia, Email nekosatriawan@gmail.com

Dear editor

We are appreciative of the authors for providing their research about “Oriented Graphene Oxide Scaffold Promotes Nerve Regeneration in vitro and in vivo” in *International Journal of Nanomedicine* 2024;19:2573–2589.¹ We would like to provide our thoughts regarding the isolation methods of dorsal root ganglia (DRG) and characterization methods of graphene oxide/polycaprolactone (GO/PCL) nanofibers as scaffolds in this study.

This study aims to create a new scaffold using graphene oxide (GO) and examine its potential therapeutic implications for nerve regeneration. In this paper, the author used trypsin to isolate DRG. Besides describing the protocol in excellent detail, the author should additionally include the references used. In earlier studies, collagenase or dispase was frequently utilized for DRG isolation in further to trypsin.²

Nanofibers are formed by overlapping layers of fibers during the electrospinning process. The number of nanofibers influences both the mass transfer from cells and the cell-substrate interaction. The thicker the layers, the lower the elastic properties of the matrix due to the density and density of the resulting nanofiber structure, which harms the infiltration of cultured cells.³

The authors did not specify the thickness of the nanofibers used in the mechanical tests and/or the length of time for electrospinning. The authors also did not explain the thickness used in mechanical tests or the thickness used for primary Schwann cells (SCs)-cultured scaffolds in vitro and in vivo. The tensile strength and Young’s modulus of the PCL nanofiber produced by electrospinning are directly correlated with its density; the higher the density, the lower the nanofiber’s Young’s modulus.⁴ The conductive nanofiber-based scaffolds are appropriate materials for neural stem cell adhesion and proliferation, but the authors did not report the conductivity value of each GO/PCL composite produced, nor did they explain the relationship between conductivity and cellular-scaffold interaction.

The author also mentioned that the combination architecture synergistically promotes nerve cell regeneration. However, the surface topography of the resulting scaffold and its influence on cell attachment are not fully described in this study. We recommend determining the porosity of the outcome nanofiber scaffold, as it can improve hydrophilicity and cell adhesion to the scaffold surface.⁵ Therefore, if using nanofiber as a scaffold, it would be required to provide a porosity value and details on the surface morphology of the resulting nanofiber. So that, it could be used as a reference for many related investigations in the future.

Acknowledgments

We should express our gratitude to Wang et al for their excellent work in this research. Furthermore, we would like to thank dr. Radiana D. Antarianto, M.Biomed., PhD., Prof. Dr. Sunarno, M.Si.Med, and all researchers at the Center for

Biomedical Research, Research Organization for Health, National Research and Innovation Agency (BRIN), for their continued assistance.

Disclosure

The authors report no conflicts of interest in this communication.

References

1. Zhou X, Tang A, Xiong C, Zhang G, Huang L, Xu F. Oriented graphene oxide scaffold promotes nerve regeneration in vitro and in vivo. *Int J Nanomed*. 2024;19:2573–2589. doi:10.2147/IJN.S439656
2. Perner C, Sokol CL. Protocol for dissection and culture of murine dorsal root ganglia neurons to study neuropeptide release. *STAR Protoc*. 2021;2(1):1–13. doi:10.1016/j.xpro.2021.100333
3. Brown JH, Das P, DiVito MD, Ivancic D, Tan LP, Wertheim JA. Nanofibrous PLGA electrospun scaffolds modified with type I collagen influence hepatocyte function and support viability in vitro. *Acta Biomater*. 2018;73:217–227. doi:10.1016/j.actbio.2018.02.009
4. Conte AA, Sun K, Hu X, Beachley VZ. Effects of fiber density and strain rate on the mechanical properties of electrospun polycaprolactone nanofiber mats. *Front Chem*. 2020;8:1–10. doi:10.3389/fchem.2020.00610
5. Lanno GM, Ramos C, Preem L, et al. Antibacterial porous electrospun fibers as skin scaffolds for wound healing applications. *ACS Omega*. 2020;5(46):30011–30022. doi:10.1021/acsomega.0c04402

Dove Medical Press encourages responsible, free and frank academic debate. The content of the International Journal of Nanomedicine 'letters to the editor' section does not necessarily represent the views of Dove Medical Press, its officers, agents, employees, related entities or the International Journal of Nanomedicine editors. While all reasonable steps have been taken to confirm the content of each letter, Dove Medical Press accepts no liability in respect of the content of any letter, nor is it responsible for the content and accuracy of any letter to the editor.

International Journal of Nanomedicine

Dovepress

Publish your work in this journal

The International Journal of Nanomedicine is an international, peer-reviewed journal focusing on the application of nanotechnology in diagnostics, therapeutics, and drug delivery systems throughout the biomedical field. This journal is indexed on PubMed Central, MedLine, CAS, SciSearch®, Current Contents®/Clinical Medicine, Journal Citation Reports/Science Edition, EMBase, Scopus and the Elsevier Bibliographic databases. The manuscript management system is completely online and includes a very quick and fair peer-review system, which is all easy to use. Visit <http://www.dovepress.com/testimonials.php> to read real quotes from published authors.

Submit your manuscript here: <https://www.dovepress.com/international-journal-of-nanomedicine-journal>

<https://doi.org/10.2147/IJN.S476121>